

## BACKGROUND

African-American (AA) women with breast cancer have a less favorable prognosis, likely due to differences in tumor biology. The Neoadjuvant Breast Cancer Symphony Trial (NBRST, NCT01479101) was a prospective trial that has shown an association of MammaPrint/Blueprint (MP/BP) with a rate of pathologic Complete Response (pCR) of 2% pCR in Luminal A with 95% Distant Metastasis Free Interval at 3 yrs. Here, we determine the MP/BP risk distribution, the response to therapy, and outcome in AA and Caucasian (Cau) patients.

## METHODS


NBRST enrolled 1,072 breast cancer pts in the US (June 2011 and December 2014). The current unplanned analysis compared clinicopathological characteristics, molecular risk assignment and outcome with neoadjuvant chemotherapy (NACT) in 150 AA and 748 Cau pts. Molecular subtyping groups were assessed by MP/BP as follows: Luminal A (MammaPrint Low Risk), Luminal B (MammaPrint High Risk), HER2 and Basal types.

## Clinical Characteristics


	AA (150)	Cau (748)	Total (898)
Med age (range)	52.5 (22-77)	54 (18-89)	<b>0.021</b>
T1-2	94 (63%)	532 (71%)	0.127
T3-4	56 (37%)	214 (29%)	
TX	0 (0%)	2 (<1%)	
N0	39 (26%)	334 (45%)	<0.001
N1	90 (60%)	331 (44%)	
N2-3	17 (11%)	55 (7%)	
NX	4 (3%)	28 (4%)	0.002
Grade 1	4 (3%)	61 (8%)	
Grade 2	40 (27%)	269 (36%)	
Grade 3	99 (66%)	397 (53%)	0.007
Grade unknown	7 (5%)	21 (3%)	
ER+	82 (55%)	496 (66%)	
PgR+	71 (47%)	406 (54%)	0.120
HER2+	40 (27%)	219 (30%)	0.476
Triple neg	48 (32%)	153 (21%)	<b>0.003</b>
NCT +/- H	145 (97%)	681 (91%)	0.059
NET +/- H	4 (3%)	61 (8%)	
NCT+ET +/- H	1 (<1%)	6 (1%)	

## MammaPrint all AA and Cau pts


### MammaPrint distribution by AA vs Cau in all pts



### MammaPrint Low Risk in all patients




### MammaPrint High Risk in all patients




## MammaPrint HR+/HER2- AA and Cau pts


### MammaPrint distribution by AA vs Cau in HR+/HER2-



### MammaPrint Low Risk in HR+/HER2- patients




### MammaPrint High Risk in HR+/HER2- patients




## Blueprint

### Blueprint distributions by AA vs Cau




### pCR rates in Blueprint Molecular Subgroups



### Uni- and multivariate analyses pCR rate

Univariate analyses	p-value	HR (95% CI)	Multivariate	p-value	HR (95% CI)
Race (AA vs Cau)	0.088	1.4 (0.95-2.0)	0.012	1.9 (1.1-3.1)	
Age	0.241		0.466		
ER	<0.001		0.797		
PR	<0.001		0.022		
HER2	<0.001		0.028		
N-stage	0.005		0.019		
T-stage	<0.001		0.001		
Grade	<0.001		<0.001		
Histol (ILC, IDC)	0.001		0.935		
BP Molecular subtype	<0.001		<0.001		

## MammaPrint in all pts




## Other tests

### TAILORx distributions and outcome

*Albain et al, SABCS 2018*  
 Distribution of risk is not different between AA and Cau in intermediate risk, although outcome of intermediate risk pts differs significantly. Risk distribution in 693 AA pts: 15% low, 68% intermediate and 17% high risk; in 8189 Cau pts: 17% low, 69% intermediate, 14% high risk. AA pts classified as intermediate risk had a 9 yr iDFS of 78.5% versus 83.9% for Cau pts (difference of 5.4%).

### PAM50 distributions and outcome

*Sun et al, SABCS 2017*  
 Racial disparities reflected by RNA based subtyping test PAM50; distribution significantly different between AA (n=454) and Cau (n=439): 31% Lum A vs 50% and 35% Basal-like vs 18%. Outcome corresponds with assigned risk:



## CONCLUSIONS

- AA pts were younger, and had a higher likelihood of having higher grade, ER-, LN+ tumors.
- AA pts were more often classified as Basal-type and less often classified as Luminal A-type compared with Cau pts by Blueprint.
- MP identified pts with Low Risk, irrespective of race (although numbers are small and longer FU is necessary).
- This study confirms that racial differences in gene expression contribute to the survival disparity observed between AA and Cau women diagnosed with breast cancer.
- MammaPrint is helpful to characterize the tumor's biology and select patients who will not benefit from chemotherapy independently of their race.