

Background

A need exists to identify biomarkers that predict chemotherapy (CT) benefit among patients with HR+HER2- early-stage breast cancer. The 70-gene signature, MammaPrint (MP), determines distant recurrence (DR) risk in early-stage breast cancer. MP has demonstrated utility in predicting adjuvant and neoadjuvant chemotherapy response¹⁻⁵. Here, we examined MP utility in predicting pathological Complete Response (pCR) to neoadjuvant CT and 5-year follow-up comparisons among patients who received adjuvant CT or endocrine therapy (ET) alone.

Methods

| Patients | Trial | Treatment | Outcome | N | |
|---|--------------------|------------------------|--------------------------------|------------|-----|
| with HR+HER2- early-stage breast cancer | NBRST | Neoadjuvant CT | pCR | 462 | |
| | FLEX (NCT03053193) | CT (+/- ET) ET only | 5-year distant recurrence risk | 181 294 | 475 |

Genomic testing

MammaPrint index is defined as UltraLow (UL: +1.000 to +0.356), Low Risk (LR: +0.355 to +0.001), High 1 (H1: 0.000 to -0.569), and High 2 (H2: -0.570 to -1.000).

Statistics

Logistic regression was used to estimate likelihood of pCR and 5-yr DR risk (including breast cancer specific death) for CT vs ET, as a continuous function of the MP index.

Results

Table 1. Clinical Characteristics of NBRST Patients

| | High 1 (N=198) | High 2 (N=129) | Low Risk (N=99) | Overall (N=426) | P-value |
|--------------------------|---------------------------|--------------------------|--------------------------|----------------------------|---------|
| Age (Years) | | | | | |
| Mean (SD) | 53 (± 12) | 50 (± 13) | 54 (± 11) | 52 (± 12) | 0.126 |
| Menopausal Status | | | | | |
| Post-Pre-/Peri- | 110 (55.6%) 88 (44.4%) | 62 (48.4%) 66 (51.6%) | 54 (55.7%) 43 (44.3%) | 226 (53.4%) 197 (46.6%) | 0.585 |
| Race | | | | | |
| Asian | 5 (2.5%) | 4 (3.1%) | 1 (1.0%) | 10 (2.3%) | 0.646 |
| Black | 27 (13.6%) | 23 (17.8%) | 7 (7.1%) | 57 (13.4%) | |
| Hispanic | 19 (9.6%) | 9 (7.0%) | 6 (6.1%) | 34 (8.0%) | |
| Other | 4 (2.0%) | 2 (1.6%) | 1 (1.0%) | 7 (1.6%) | |
| White | 143 (72.2%) | 91 (70.5%) | 84 (84.8%) | 318 (74.6%) | |
| Tumor Stage | | | | | |
| T1 | 22 (11.2%) | 12 (9.4%) | 12 (12.1%) | 46 (10.9%) | 0.964 |
| T2 | 108 (54.8%) | 75 (59.1%) | 51 (51.5%) | 234 (55.3%) | |
| T3 | 54 (27.4%) | 30 (23.6%) | 31 (31.3%) | 115 (27.2%) | |
| T4 | 13 (6.6%) | 10 (7.9%) | 5 (5.1%) | 28 (6.6%) | |
| Lymph Node Stage | | | | | |
| N0 | 56 (29.6%) | 43 (35.2%) | 39 (41.1%) | 138 (34.0%) | 0.476 |
| N1 | 115 (60.8%) | 61 (50.0%) | 46 (48.4%) | 222 (54.7%) | |
| N2 | 16 (8.5%) | 14 (11.5%) | 7 (7.4%) | 37 (9.1%) | |
| N3 | 2 (1.1%) | 4 (3.3%) | 3 (3.2%) | 9 (2.2%) | |
| Grade | | | | | |
| G1 | 15 (7.9%) | 2 (1.6%) | 22 (23.7%) | 39 (9.6%) | <0.001 |
| G2 | 95 (50.0%) | 18 (14.6%) | 57 (61.3%) | 170 (41.9%) | |
| G3 | 80 (42.1%) | 103 (83.7%) | 14 (15.1%) | 197 (48.5%) | |

Data presented in both tables as N (%) unless indicated otherwise; Unknown values excluded; N = sample size; SD = standard deviation

- Neoadjuvantly treated patients were significantly more likely to achieve pCR as MP risk increased, with up to 30% pCR observed in H2 tumors ($p < 0.001$) (**Figure 1**).
- The MP index demonstrated effective performance for predicting 5-yr DR risk in the ET group, with an Area under the ROC Curve (AUC) of 0.73 ($p=0.005$), and the CT group, with an AUC of 0.77 ($p=0.008$) (**Figure 2**).
- ET only treated pts had greater 5-yr DR risk with increasing MP risk compared to patients treated with CT.
- Comparing ET vs CT groups, <1.5% and <1.0% difference in 5-yr DR risk was observed for MP indices in the LR and UL range, respectively.
- In contrast, CT benefit increased with increasing MP risk, with 2-7% absolute risk difference observed in H1, and 8-16% difference among H2 tumors.

Figure 1. Probability of achieving pCR with Neoadjuvant Chemotherapy across MammaPrint Index





Table 2. Clinical Characteristics of FLEX Patients

| | CT +/- ET (N=181) | ET only (N=294) | Overall (N=475) | P-value |
|----------------------------------|---------------------------|---------------------------|---------------------------|---------|
| Age (Years) | | | | |
| Mean (SD) | 59 (± 11) | 62 (± 11) | 61 (± 11) | 0.0597 |
| Menopausal Status | | | | |
| Post-Pre-/Peri- | 134 (75.7%) 43 (24.3%) | 239 (84.5%) 44 (15.5%) | 373 (81.1%) 87 (18.9%) | 0.0735 |
| Race | | | | |
| Asian or Pacific Islander | 4 (2.3%) | 5 (1.8%) | 9 (2.0%) | 0.252 |
| American Indian or Alaska Native | 0 (0%) | 1 (0.4%) | 1 (0.2%) | |
| Black | 12 (6.8%) | 12 (4.3%) | 24 (5.3%) | |
| Latin American/Hispanic | 0 (0%) | 10 (3.6%) | 10 (2.2%) | |
| White | 160 (90.9%) | 253 (90.0%) | 413 (90.4%) | |
| Tumor Stage | | | | |
| T1 | 86 (64.7%) | 124 (68.9%) | 210 (67.1%) | 0.366 |
| T2 | 39 (29.3%) | 53 (29.4%) | 92 (29.4%) | |
| T3 | 8 (6.0%) | 3 (1.7%) | 11 (3.5%) | |
| Lymph Node Status | | | | |
| Negative | 104 (83.2%) | 153 (87.9%) | 257 (86.0%) | 0.526 |
| Positive | 21 (16.8%) | 21 (12.1%) | 42 (14.0%) | |
| Lymph Node Stage | | | | |
| N0 | 104 (83.2%) | 153 (87.9%) | 257 (86.0%) | 0.769 |
| N1 | 20 (16.0%) | 19 (10.9%) | 39 (13.0%) | |
| N2 | 1 (0.8%) | 2 (1.1%) | 3 (1.0%) | |
| N3 | 0 (0%) | 0 (0%) | 0 (0%) | |
| Grade | | | | |
| G1 | 46 (26.4%) | 91 (32.3%) | 137 (30.0%) | 0.777 |
| G2 | 98 (56.3%) | 147 (52.1%) | 245 (53.7%) | |
| G3 | 30 (17.2%) | 44 (15.6%) | 74 (16.2%) | |

Figure 2. Probability of Distant Recurrence with Adjuvant CT or ET alone across MammaPrint Index



Conclusions

- These data show that MammaPrint indices within LR and UL ranges exhibit low chemosensitivity and do not derive significant CT benefit, consistent with results in MINDACT¹.
- The increasing chemotherapy benefit observed with increasing MP risk is consistent with the reduction in recurrence risk reported by Knauer et al.² for adjuvant therapy, as well as neoadjuvant results from I-SPY 2, NBRST, and FLEX data³⁻⁵.
- These data indicate the utility of MammaPrint to predict neoadjuvant and adjuvant CT benefit in patients with HR+HER2- breast cancer.